



Wireless Elapse Timers



This is a typical CB625-RS digital clock with an elapse timer.



BRG Elapse Timers are designed to operate in the sterile conditions of an operating room. The timers can either be stand alone or can function as a timer and a digital clock. This gives you a synchronized time piece in the operating room and a timer at the touch of a button. The timers are designed to either count up or countdown.

The timers can be controlled by one of five methods:

1. Control buttons mounted on the timer.
2. Infrared remote control.
3. Wired remote control.
4. Wall mounted controls.
5. Wired connection box to connect your control switches.

The LED character size can be 1.8", 2.5" or 4.0" in height.



This is a typical CB425-RS with elapse timer



600 N River Street • Derby, Kansas 67037
800-295-0220 • 316-788-2000
Sales Fax: 720-293-9393
Email: sales@brgproducts.com
Web Site: <http://www.brgwireless.com>